

EBP12320 - IMMUNIZING PEPTIDE

SKU: EBP12320



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001209.1; NP_996808.1; NP_001280571.1
Immunogen	Peptide with sequence C-SHLQHVKYKGQE, from the internal region of the protein sequence according to NP_001209.1; NP_996808.1; NP_001280571.1.
Product	The peptide was used in the production of EBP12320 - Goat Anti-carbonic anhydrase XII (aa188-199) Antibody.
Comments	It is currently untested in blocking applications.
Peptide Sequence	C-SHLQHVKYKGQE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)