

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12326 - IMMUNIZING PEPTIDE

SKU: EBP12326

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_000083.3 ID

Peptide with sequence C-SEEAIRPYVSR, from the internal region (near C terminus) of the protein sequence

according to NP_000083.3.

The peptide was used in the production of EB12326 - Goat Anti-COL4A4 (aa1555-1565) Antibody. It is currently Product

Comments untested in blocking applications.

Peptide

C-SEEAIRPYVSR Sequence

Shipping

Ambient Instructions

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png Grade