

EBP12326 - IMMUNIZING PEPTIDE

SKU: EBP12326



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_000083.3

Immunogen Peptide with sequence C-SEEAIRPYVSR, from the internal region (near C terminus) of the protein sequence according to NP_000083.3.

Product Comments The peptide was used in the production of EB12326 - Goat Anti-COL4A4 (aa1555-1565) Antibody. It is currently untested in blocking applications.

Peptide Sequence C-SEEAIRPYVSR

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)