

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP12340 - IMMUNIZING PEPTIDE**

**SKU:** EBP12340

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage Shipped at ambient temperature, store at -20°C.

**Accession** 

NP\_004235.4; NP\_981961.2; ID

Peptide with sequence C-SRGENLVHQIQYRE, from the internal region (near C terminus) of the protein sequence

according to NP\_004235.4; NP\_981961.2;.

The peptide was used in the production of EB12340 - Goat Anti-CD163 Antibody. It is currently untested in blocking **Product** 

Comments applications.

**Peptide** 

C-SRGENLVHQIQYRE Sequence

Shipping

Ambient Instructions

**Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite\_plus\_medium.png$ Grade