

EBP12340 - IMMUNIZING PEPTIDE

SKU: EBP12340



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_004235.4; NP_981961.2;
Immunogen	Peptide with sequence C-SRGENLVHQIQYRE, from the internal region (near C terminus) of the protein sequence according to NP_004235.4; NP_981961.2;.
Product Comments	The peptide was used in the production of EB12340 - Goat Anti-CD163 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-SRGENLVHQIQYRE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png