

EBP12372 - IMMUNIZING PEPTIDE

SKU: EBP12372



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_058197.2; NP_002526.2; NP_001027903.1
Immunogen	Peptide with sequence C-RKTVLRGNSDGT, from the internal region (near N terminus) of the protein sequence according to NP_058197.2; NP_002526.2; NP_001027903.1.
Product	The peptide was used in the production of EB12372 - Goat Anti-OAS2 Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	C-RKTVLRGNSDGT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)