

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12376 - IMMUNIZING PEPTIDE

SKU: EBP12376

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 000499.1; NP 068657.1 ID

Peptide with sequence C-STSYNRGDSTFES, from the internal region (near C terminus) of the protein sequence **Immunogen**

according to NP 000499.1; NP 068657.1.

The peptide was used in the production of EB12376 - Goat Anti-fibrinogen alpha chain Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-STSYNRGDSTFES Sequence

Shipping

Grade

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png