

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12378 - IMMUNIZING PEPTIDE

SKU: EBP12378

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

ID

NP 005132.2; NP 001171670.1

Peptide with sequence C-EKHQLYIDETVNSN, from the internal region of the protein sequence according to Immunogen

NP 005132.2; NP 001171670.1.

The peptide was used in the production of EB12378 - Goat Anti-fibrinogen beta chain (aa177-190) Antibody. It **Product**

Comments is currently untested in blocking applications.

Peptide

C-EKHQLYIDETVNSN Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade