

EBP12380 - IMMUNIZING PEPTIDE

SKU: EBP12380



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_000500.2; NP_068656.2
Immunogen	Peptide with sequence C-KDTVQIHDITGKD, from the internal region of the protein sequence according to NP_000500.2; NP_068656.2.
Product	The peptide was used in the production of EB12380 - Goat Anti-fibrinogen gamma chain (aa166-178) Antibody.
Comments	It is currently untested in blocking applications.
Peptide Sequence	C-KDTVQIHDITGKD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)