

EBP12386 - IMMUNIZING PEPTIDE

SKU: EBP12386



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_079409.3; NP_001177203.1; XP_005264629.1; XP_005264630.1; XP_005264631.1; XP_005264632.1
Immunogen	Peptide with sequence C-KISSYTSYPMHD, from the internal region of the protein sequence according to NP_079409.3; NP_001177203.1; XP_005264629.1; XP_005264630.1; XP_005264631.1; XP_005264632.1.
Product Comments	The peptide was used in the production of EB12386 - Goat Anti-FBXO11 (aa754-765) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KISSYTSYPMHD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)