

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP12387 - IMMUNIZING PEPTIDE**

**SKU:** EBP12387

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

**Accession** 

ID

NP\_004334.1

 $\begin{tabular}{l|l} \textbf{Immunogen} \\ \hline & \textbf{Peptide with sequence C-HPEIDNPEYSPDPS, from the internal region of the protein sequence according to NP\_004334.1.} \\ \end{tabular}$ 

**Product** The peptide was used in the production of EB12387 - Goat Anti-calreticulin Antibody. It is currently untested in

**Comments** blocking applications.

**Peptide** 

Sequence

C-HPEIDNPEYSPDPS

**Shipping** 

Ambient Instructions

**Product** 

Grade

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$