

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP12389 - IMMUNIZING PEPTIDE**

**SKU:** EBP12389

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 004037.1; NP 001244263.1; NP 001001935.1 ID

Peptide with sequence C-RVHGLRNVQAEE, from the internal region (near N terminus) of the protein sequence **Immunogen** 

according to NP 004037.1; NP 001244263.1; NP 001001935.1.

The peptide was used in the production of EB12389 - Goat Anti-ATP5A1 Antibody. It is currently untested in **Product** 

Comments blocking applications.

**Peptide** C-RVHGLRNVQAEE Sequence

Shipping

Instructions Ambient

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade