

# EBP12392 - IMMUNIZING PEPTIDE

**SKU:** EBP12392



---

## SPECIFICATIONS

<b>Formulation</b>	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
<b>Unit Size</b>	100 µg
<b>Storage</b>	Shipped at ambient temperature, store at -20°C.
<b>Instructions</b>	
<b>Accession ID</b>	NP_000025.1; NP_001230106.1
<b>Immunogen</b>	Peptide with sequence C-QKADDGRFPFQ, from the internal region of the protein sequence according to NP_000025.1; NP_001230106.1.
<b>Product Comments</b>	The peptide was used in the production of EB12392 - Goat Anti-ALDOA (aa86-96) Antibody. It is currently untested in blocking applications.
<b>Peptide Sequence</b>	C-QKADDGRFPFQ
<b>Shipping Instructions</b>	Ambient
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>

## DOCUMENTS

- [Data Sheet](#)