

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12411 - IMMUNIZING PEPTIDE

SKU: EBP12411

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 004823.1; NP 001177931.1; NP 001177932.1 ID

Peptide with sequence C-KEDPTVSALLTSEKD, from the internal region (near C terminus) of the protein **Immunogen**

sequence according to NP 004823.1; NP 001177931.1; NP 001177932.1.

The peptide was used in the production of EB12411 - Goat Anti-GSTO1 Antibody. It is currently untested in **Product**

Comments blocking applications. **Peptide**

C-KEDPTVSALLTSEKD Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade