

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP12417 - IMMUNIZING PEPTIDE**

**SKU:** EBP12417

## **SPECIFICATIONS**

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$ 

**Unit Size** 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

**Accession** 

NP 001138632.1; NP 003281.1 ID

Peptide with sequence C-ENVGLHQTLDQ, from the C Terminus of the protein sequence according to Immunogen

NP 001138632.1; NP 003281.1.

The peptide was used in the production of EB12417 - Goat Anti-tropomyosin 4 / TPM4 Antibody. It is currently **Product** 

**Comments** untested in blocking applications.

**Peptide** 

C-ENVGLHQTLDQ Sequence

Shipping

Instructions Ambient

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png