

EBP12420 - IMMUNIZING PEPTIDE

SKU: EBP12420



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001123476.1; NP_001093.1; NP_001123477.1
Immunogen	Peptide with sequence C-TPQEINGKWDHVRQ, from the internal region of the protein sequence according to NP_001123476.1; NP_001093.1; NP_001123477.1.
Product	The peptide was used in the production of EB12420 - Goat Anti-ACTN1 (aa596-609) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-TPQEINGKWDHVRQ
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)