

EBP12421 - IMMUNIZING PEPTIDE

SKU: EBP12421



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001738.2; NP_001243069.1; NP_001307663.1
Immunogen	Peptide with sequence C-NIRATERALNWDS, from the internal region of the protein sequence according to NP_001738.2; NP_001243069.1; NP_001307663.1.
Product	The peptide was used in the production of EB12421 - Goat Anti-CAPG Antibody. It is currently untested in
Comments	blocking applications.
Peptide Sequence	C-NIRATERALNWDS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)