

EBP12445 - IMMUNIZING PEPTIDE

SKU: EBP12445



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_443170.1; NP_001152870.1
Immunogen	Peptide with sequence C-SRLRKANITDVD, from the C Terminus of the protein sequence according to NP_443170.1; NP_001152870.1.
Product Comments	The peptide was used in the production of EBP12445 - Goat Anti-FCRL1 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-SRLRKANITDVD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)