

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12487 - IMMUNIZING PEPTIDE

SKU: EBP12487

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

Accession

NP_663610.1 ID

Immunogen Peptide with sequence CEETFPGEEDGLH, from the C Terminus of the protein sequence according to NP_663610.1.

Product The peptide was used in the production of EB12487 - Goat Anti-Adipogenin (mouse) Antibody. It is currently

Comments untested in blocking applications.

Peptide

CEETFPGEEDGLH Sequence

Shipping

Ambient Instructions **Product**

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png