

EBP12488 - IMMUNIZING PEPTIDE

SKU: EBP12488



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_663610.1

Immunogen Peptide with sequence SQESKEDDSDLC, from the internal region of the protein sequence according to NP_663610.1.

Product Comments The peptide was used in the production of EBP12488 - Goat Anti-Adipogenin (mouse, aa41-52) Antibody. It is currently untested in blocking applications.

Peptide Sequence SQESKEDDSDLC

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)