

EBP12501 - IMMUNIZING PEPTIDE

SKU: EBP12501



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001768.1; NP_942088.1; NP_001369235.1
Immunogen	Peptide with sequence C-SNQKTIQPPRK, from the internal region (near C terminus) of the protein sequence according to NP_001768.1; NP_942088.1; NP_001369235.1.
Product Comments	The peptide was used in the production of EB12501 - Goat Anti-CD47 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-SNQKTIQPPRK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)