

EBP12515 - IMMUNIZING PEPTIDE

SKU: EBP12515



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001164645.1; NP_001328.1
Immunogen	Peptide with sequence C-DQFPESVTENFEYD, from the internal region (near N terminus) of the protein sequence according to NP_001164645.1; NP_001328.1.
Product Comments	The peptide was used in the production of EBP12515 - Goat Anti-CX3CR1 (aa34-47) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-DQFPESVTENFEYD
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)