

EBP12527 - IMMUNIZING PEPTIDE

SKU: EBP12527



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_009168.4; NP_001005176.1; NP_001265380.1; NP_001265381.1; NP_001265382.1
Immunogen	Peptide with sequence EGQNLQEQVCPEP, from the N-Terminus of the protein sequence according to NP_009168.4; NP_001005176.1; NP_001265380.1; NP_001265381.1; NP_001265382.1.
Product Comments	The peptide was used in the production of EB12527 - Goat Anti-SP140 (aa29-41) Antibody. It is currently untested in blocking applications.
Peptide Sequence	EGQNLQEQVCPEP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Safety Data Sheet](#)