

EBP12531 - IMMUNIZING PEPTIDE

SKU: EBP12531



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_620412.1; NP_733839.1
Immunogen	Peptide with sequence C-SPDSLRYFSK, from the internal region of the protein sequence according to NP_620412.1; NP_733839.1.
Product Comments	The peptide was used in the production of EBP12531 - Goat Anti-MSI2 / musashi-2 (aa33-43) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-SPDSLRYFSK
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)