

EBP12563 - IMMUNIZING PEPTIDE

SKU: EBP12563



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_055181.3; NP_996561.1
Immunogen	Peptide with sequence C-TENKAKSGKQNS, from the internal region of the protein sequence according to NP_055181.3; NP_996561.1.
Product Comments	The peptide was used in the production of EBP12563 - Goat Anti-GNL3 (aa115-126) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-TENKAKSGKQNS
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)