



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12625 - IMMUNIZING PEPTIDE

SKU: EBP12625

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage

Shipped at ambient temperature, store at -20°C.

Accession

NP_065434.1; NP_001158162.1; NP_001158163.1 ID

Immunogen Peptide with sequence ESLHCDPSPRAH, from the internal region of the protein sequence according to NP_065434.1; NP_001158162.1; NP_001158163.1.

The peptide was used in the production of EB12625 - Goat Anti-AMHR2 (aa88-99) Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

ESLHCDPSPRAH Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws. com/assets/products/image/elite_plus_medium.png$ Grade