

EBP12639 - IMMUNIZING PEPTIDE

SKU: EBP12639



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_002650.1; NP_001005376.1; NP_001005377.1
Immunogen	Peptide with sequence C-RTGLKITSLTE , from the internal region of the protein sequence according to NP_002650.1; NP_001005376.1; NP_001005377.1.
Product	The peptide was used in the production of EBP12639 - Goat Anti-U-PAR / CD87 Antibody. It is currently untested
Comments	in blocking applications.
Peptide Sequence	C-RTGLKITSLTE
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)