

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12639 - IMMUNIZING PEPTIDE

SKU: EBP12639

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

ID

NP 002650.1; NP 001005376.1; NP 001005377.1

Peptide with sequence C-RTGLKITSLTE , from the internal region of the protein sequence according to Immunogen

NP 002650.1; NP 001005376.1; NP 001005377.1.

The peptide was used in the production of EB12639 - Goat Anti-U-PAR / CD87 Antibody. It is currently untested **Product**

Comments in blocking applications.

Peptide

C-RTGLKITSLTE Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade