

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **EBP12715 - IMMUNIZING PEPTIDE**

**SKU:** EBP12715

## **SPECIFICATIONS**

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

**Unit Size** 100 μg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

NP 066978.1; NP 001138368.1; NP 001138369.1; NP 001138365.1; NP 001138366.1; NP 001138367.1; Accession

ID NP 001138371.1

Peptide with sequence C-NRGEPNNVGEED, from the internal region of the protein sequence according to

Immunogen NP\_066978.1; NP\_001138368.1; NP\_001138369.1; NP\_001138365.1; NP\_001138366.1; NP\_001138367.1;

NP\_001138371.1.

Product The peptide was used in the production of EB12715 - Goat Anti-DC-SIGN / CD209 Antibody. It is currently untested

Comments in blocking applications.

**Peptide** 

Sequence

C-NRGEPNNVGEED

Shipping

Ambient Instructions

**Product** Grade

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_plus\_medium.png$