

EBP12749 - IMMUNIZING PEPTIDE

SKU: EBP12749



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001124459.1; NP_001123927.1; NP_001124458.1; NP_001124457.1; NP_001124456.1; NP_001124455.1; NP_001124454.1; NP_003485.1; NP_001124448.1; NP_001124449.1; NP_001124450.1; NP_001124451.1; NP_001124452.1; NP_001124453.1
Immunogen	Peptide with sequence HLFCQQHRVKAP, from the internal region of the protein sequence according to NP_001124459.1; NP_001123927.1; NP_001124458.1; NP_001124457.1; NP_001124456.1; NP_001124455.1; NP_001124454.1; NP_003485.1; NP_001124448.1; NP_001124449.1; NP_001124450.1; NP_001124451.1; NP_001124452.1; NP_001124453.1.
Product Comments	The peptide was used in the production of EBP12749 - Goat Anti-dysferlin Antibody. It is currently untested in blocking applications.
Peptide Sequence	HLFCQQHRVKAP
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

DOCUMENTS

- [Data Sheet](#)