



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12749 - IMMUNIZING PEPTIDE

SKU: EBP12749

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage

Instructions Shipped at ambient temperature, store at -20°C.

NP 001124459.1; NP 001123927.1; NP 001124458.1; NP 001124457.1; NP 001124456.1; NP 001124455.1;

Accession NP 001124454.1; NP 003485.1; NP 001124448.1; NP 001124449.1; NP 001124450.1; NP 001124451.1;

ID NP_001124452.1; NP_001124453.1

Peptide with sequence HLFCQQHRVKAP, from the internal region of the protein sequence according to NP_001124459.1; NP_001123927.1; NP_001124458.1; NP_001124457.1; NP_001124456.1; NP_001124455.1;

Immunogen NP_001124454.1; NP_003485.1; NP_001124448.1; NP_001124449.1; NP_001124450.1; NP_001124451.1;

NP_001124452.1; NP_001124453.1.

The peptide was used in the production of EB12749 - Goat Anti-dysferlin Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

HLFCQQHRVKAP Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade