

EBP12762 - IMMUNIZING PEPTIDE

SKU: EBP12762



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_001970.2; NP_001243411.1; NP_001243412.1
Immunogen	Peptide with sequence C-SRTFKSLFRASDES, from the internal region of the protein sequence according to NP_001970.2; NP_001243411.1; NP_001243412.1.
Product Comments	The peptide was used in the production of EBP12762 - Goat Anti-EPHX2 Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-SRTFKSLFRASDES
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)