

EBP12801 - IMMUNIZING PEPTIDE

SKU: EBP12801



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NP_001499.1

Immunogen Peptide with sequence C-HAHSTTDSARFVQRP, from the internal region of the protein sequence according to NP_001499.1.

Product Comments The peptide was used in the production of EBP12801 - Goat Anti-GPR39 (aa377-391) Antibody. It is currently untested in blocking applications.

Peptide Sequence C-HAHSTTDSARFVQRP

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)