

EBP12830 - IMMUNIZING PEPTIDE

SKU: EBP12830



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_000181.2; NP_001019553.1; NP_001245137.1; NP_001245138.1
Immunogen	Peptide with sequence C-HLEFRSIRGNLNT, from the internal region of the protein sequence according to NP_000181.2; NP_001019553.1; NP_001245137.1; NP_001245138.1.
Product	The peptide was used in the production of EBP12830 - Goat Anti-HMBS (aa160-173) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-HLEFRSIRGNLNT
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)