

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12834 - IMMUNIZING PEPTIDE

SKU: EBP12834

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession

NP_002147.2 ID

 $\begin{tabular}{ll} \textbf{Immunogen} \\ \textbf{NP_002147.2}. \end{tabular} Peptide with sequence C-NLEDVQPHDLGK, from the internal region of the protein sequence according to NP_002147.2. \end{tabular}$

The peptide was used in the production of EB12834 - Goat Anti-HSP60 / HSPD1 (aa333-344) Antibody. It is currently **Product**

Comments untested in blocking applications.

Peptide

C-NLEDVQPHDLGK Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade