

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP12856 - IMMUNIZING PEPTIDE

SKU: EBP12856

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

100 μg **Unit Size**

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NP 005338.1 ID

Peptide with sequence C-KVVEKKTKPYIQVD, from the internal region of the protein sequence according to Immunogen

NP 005338.1.

The peptide was used in the production of EB12856 - Goat Anti-HSPA5 / GRP78 (aa118-131) Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-KVVEKKTKPYIQVD Sequence

Shipping

Instructions Ambient

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade