

EBP12966 - IMMUNIZING PEPTIDE

SKU: EBP12966



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_803228.2; NP_663482.2
Immunogen	Peptide with sequence SRKGQEKPREIVD-C, from the N Terminus of the protein sequence according to NP_803228.2; NP_663482.2.
Product	The peptide was used in the production of EBP12966 - Goat Anti-Aspartyl-tRNA Synthetase (N terminus)
Comments	Antibody. It is currently untested in blocking applications.
Peptide Sequence	SRKGQEKPREIVD-C
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)