

EBP12988 - IMMUNIZING PEPTIDE

SKU: EBP12988



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_003670.2
Immunogen	Peptide with sequence C-KRNFQIDVYEARED, from the internal region (near N terminus) of the protein sequence according to NP_003670.2.
Product Comments	The peptide was used in the production of EBP12988 - Goat Anti-kynurenine 3-monooxygenase (aa29-42) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-KRNFQIDVYEARED
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)