

EBP12997 - IMMUNIZING PEPTIDE

SKU: EBP12997



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_001026971.1; NP_005560.1; NP_057952.1
Immunogen	Peptide with sequence C-QVKEVNRMHISPN, from the internal region of the protein sequence according to NP_001026971.1; NP_005560.1; NP_057952.1.
Product Comments	The peptide was used in the production of EB12997 - Goat Anti-LIMK2 (aa181-193) Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-QVKEVNRMHISPN
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

DOCUMENTS

- [Data Sheet](#)