

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP13088 - IMMUNIZING PEPTIDE

SKU: EBP13088

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

PODTC2; NC 045512.2/MN908947 ID

Peptide with sequence PROPRIETARY-SARS-CoV-2, from the N Terminus of the protein sequence according to Immunogen

P0DTC2; NC_045512.2/MN908947.

The peptide was used in the production of EB13088 - Goat Anti-SARS-COV2 Spike protein- S1 Domain (N **Product**

Comments Terminal) Antibody. It is currently untested in blocking applications.

Peptide

PROPRIETARY-SARS-CoV-2 Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png