

EBP13088 - IMMUNIZING PEPTIDE

SKU: EBP13088



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID P0DTC2; NC_045512.2/MN908947

Immunogen Peptide with sequence PROPRIETARY-SARS-CoV-2, from the N Terminus of the protein sequence according to P0DTC2; NC_045512.2/MN908947.

Product Comments The peptide was used in the production of EBP13088 - Goat Anti-SARS-COV2 Spike protein- S1 Domain (N Terminal) Antibody. It is currently untested in blocking applications.

Peptide Sequence PROPRIETARY-SARS-CoV-2

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png