

EBP13089 - IMMUNIZING PEPTIDE

SKU: EBP13089



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NC_045512.2/MN908947
Immunogen	Peptide with sequence C-PROPRIETARY.., from the C or N Terminus of the protein sequence according to NC_045512.2/MN908947.
Product	The peptide was used in the production of EB13089 - Goat Anti-Spike protein- S1 Domain (aa672-682)
Comments	Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PROPRIETARY..
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png