

EBP13090 - IMMUNIZING PEPTIDE

SKU: EBP13090



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	
Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NC_045512.2/MN908947
Immunogen	Peptide with sequence C-PROPRIETARY,,, from the C or N Terminus of the protein sequence according to NC_045512.2/MN908947.
Product	The peptide was used in the production of EBP13090 - Goat Anti-Spike protein- S2 Domain (aa781-791)
Comments	Antibody. It is currently untested in blocking applications.
Peptide Sequence	C-PROPRIETARY,,
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png