

EBP13091 - IMMUNIZING PEPTIDE

SKU: EBP13091



SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 µg

Storage Instructions Shipped at ambient temperature, store at -20°C.

Accession ID NC_045512.2/MN908947

Immunogen Peptide with sequence C-PROPRIETARY,, from the C or N Terminus of the protein sequence according to NC_045512.2/MN908947.

Product The peptide was used in the production of EBP13091 - Goat Anti-Spike protein- S2 Domain (aa810-820)

Comments Antibody. It is currently untested in blocking applications.

Peptide Sequence C-PROPRIETARY,.

Shipping Instructions Ambient

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png