

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP13100 - IMMUNIZING PEPTIDE

SKU: EBP13100

SPECIFICATIONS

 $\textbf{Formulation} \ 100 \mu g \ of \ Lyophilised \ peptide, \ formulation \ De-ionized \ (DI) \ water.$

Unit Size 100 μα

Storage

Shipped at ambient temperature, store at -20°C. Instructions

Accession

NC 045512.2; MN908947 ID

Peptide with sequence C-10 aa long, RBD of spike protein, from the internal region of the protein sequence **Immunogen**

according to NC 045512.2; MN908947.

The peptide was used in the production of EB13100 - Goat Anti-SARS CoV2 Spike protein-RBD Antibody. It is **Product**

Comments currently untested in blocking applications.

Peptide

C-10 aa long, RBD of spike protein Sequence

Shipping

Instructions Ambient

Product Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png