

EBP13100 - IMMUNIZING PEPTIDE

SKU: EBP13100



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NC_045512.2; MN908947
Immunogen	Peptide with sequence C-10 aa long, RBD of spike protein, from the internal region of the protein sequence according to NC_045512.2; MN908947.
Product	The peptide was used in the production of EBP13100 - Goat Anti-SARS CoV2 Spike protein-RBD Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	C-10 aa long, RBD of spike protein
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png