

EBP13106 - IMMUNIZING PEPTIDE

SKU: EBP13106



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage Instructions	Shipped at ambient temperature, store at -20°C.
Accession ID	NP_037393.1, NP_001317680.1, NP_001317681.1, NP_001317682.1 , NP_001341756.1
Immunogen	Peptide with sequence RDSVSPPKSLFSQC, from the internal region of the protein sequence according to NP_037393.1, NP_001317680.1, NP_001317681.1, NP_001317682.1 , NP_001341756.1.
Product Comments	The peptide was used in the production of EBP13106 - Goat Anti-PPARGC1A Antibody. It is currently untested in blocking applications.
Peptide Sequence	RDSVSPPKSLFSQC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png