

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

EBP13106 - IMMUNIZING PEPTIDE

SKU: EBP13106

SPECIFICATIONS

Formulation 100µg of Lyophilised peptide, formulation De-ionized (DI) water.

Unit Size 100 μg

Storage Shipped at ambient temperature, store at -20°C.

Accession $NP_037393.1,\,NP_001317680.1,\,NP_001317681.1,\,NP_001317682.1\,,\,NP_001341756.1$ ID

 $\textbf{Immunogen} \begin{tabular}{l} \textbf{Peptide with sequence RDSVSPPKSLFSQC}, from the internal region of the protein sequence according to NP_037393.1, NP_001317680.1, NP_001317681.1, NP_001317682.1, NP_001341756.1. \end{tabular}$

The peptide was used in the production of EB13106 - Goat Anti-PPARGC1A Antibody. It is currently untested in **Product**

Comments blocking applications.

Peptide

RDSVSPPKSLFSQC Sequence

Shipping

Ambient Instructions

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_plus_medium.png$ Grade