

EBP13107 - IMMUNIZING PEPTIDE

SKU: EBP13107



SPECIFICATIONS

Formulation	100µg of Lyophilised peptide, formulation De-ionized (DI) water.
Unit Size	100 µg
Storage	Shipped at ambient temperature, store at -20°C.
Instructions	
Accession ID	NP_000659.2; NP_000658.1; NP_000660.1
Immunogen	Peptide with sequence STAGKVIKCKAC, from the N Terminus of the protein sequence according to NP_000659.2; NP_000658.1; NP_000660.1.
Product	The peptide was used in the production of EBP13107 - Goat Anti-pan ADH (N terminus) Antibody. It is currently untested in blocking applications.
Comments	
Peptide Sequence	STAGKVIKCKAC
Shipping Instructions	Ambient
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png