

GOAT ANTI-14-3-3 TAU / YWHAQ ANTIBODY

SKU: EB08137

250kDa
150kDa
100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	tyrosine 3/tryptophan 5 -monooxygenase activation protein, theta polypeptide protein tau 14-3-3 protein theta 14-3-3 protein T-cell HS1 1C5 14-3-3 tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, theta polypeptide 14-3-3 protein tau YWHAQ
Accession ID	NP_006817.1
Blocking Peptide	EBP08137
Immunogen	Peptide with sequence C-DDRKQTIDNSQ, from the internal region of the protein sequence according to NP_006817.1.
Peptide Sequence	C-DDRKQTIDNSQ
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	10971
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:4000.
Western Blot	Approx 28kDa band observed in Human Brain (Hippocampus) lysates (calculated MW of 27.8kDa according to NP_006817.1). Recommended concentration: 0.03-0.1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

